


E 109

Insert

for the optimum loading of beakers up to 250 ml.



- Suitable for 21 glass beakers, max. 250 mL
- equipped with 21 x 3 spikes
- For use in upper and lower basket

EAN: 4002511115729 / Material number: 03808360 / Old Material Number: 69510901D

Product category	
Glassware Washers, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 6011	•
PLW 6111	•
PLW 8505	•
Product category	
Insert / module for beakers	•
Product characteristics	
Material	Stainless steel, Wire, rilsan-coated
Colour	Stainless Steel, White
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Beakers [quantity]	21
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	155
External dimensions, net width in mm	230
External dimensions, net depth in mm	460
External dimensions, gross height in mm	7 5/16 (185)
External dimensions, gross width in mm	9 1/16 (230)
External dimensions, gross depth in in (mm)	19 5/16 (490)
Net weight in kg	3 (1.21)
Gross weight in kg	4 (1.54)