


E 124

Insert

For the optimum loading of infusion flasks. Surface area 100 x 100 mm.



- 16 compartments, dimension: 100 x 100 mm (bottle), 48 x 48 mm (neck)
- For use in lower basket

EAN: 4002511115811 / Material number: 03808500 / Old Material Number: 69512401D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Insert for various bottles	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Number of compartments	16
Application area	
Reprocessing of surgical instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	148
External dimensions, net width in mm	460
External dimensions, net depth in mm	460
External dimensions, gross height in mm	6 3/4 (170)
External dimensions, gross width in mm	18 15/16 (480)
External dimensions, gross depth in in (mm)	18 15/16 (480)
Net weight in kg	4 (1.44)
Gross weight in kg	4 (1.80)