


E 125

Insert

For the optimum loading of infusion flasks. Surface area 125 x 125 mm.



- 9 compartments, dimension: 125 x 125 mm (bottle), 55 x 55 mm (neck)
- For use in lower basket

EAN: 4002511115828 / Material number: 03808510 / Old Material Number: 69512501D

<b>Product category</b>	
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
<b>Product category</b>	
Insert for various bottles	•
<b>Product characteristics</b>	
Material	Stainless steel
Colour	Stainless Steel
Number of compartments	9
<b>Application area</b>	
Reprocessing of surgical instruments	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	224
External dimensions, net width in mm	460
External dimensions, net depth in mm	460
External dimensions, gross height in mm	9 11/16 (245)
External dimensions, gross width in mm	18 15/16 (480)
External dimensions, gross depth in in (mm)	18 15/16 (480)
Net weight in kg	5 (1.81)
Gross weight in kg	6 (2.56)