


E 126

Insert

For the optimum loading of infusion flasks. Surface area 45 x 45 mm.



- 48 compartments, dimension: 45 x 45 mm (bottle), 28 x 28 mm (neck)

EAN: 4002511115835 / Material number: 03808520 / Old Material Number: 69512601D

<b>Product category</b>	
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
<b>Product category</b>	
Insert for various bottles	•
<b>Product characteristics</b>	
Material	Stainless steel
Colour	Stainless Steel
Number of compartments	48
<b>Application area</b>	
Reprocessing of surgical instruments	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	83
External dimensions, net width in mm	445
External dimensions, net depth in mm	445
External dimensions, gross height in mm	6 5/16 (160)
External dimensions, gross width in mm	18 3/4 (475)
External dimensions, gross depth in in (mm)	18 3/4 (475)
Net weight in kg	5 (1.92)
Gross weight in kg	5 (2.05)