

E 129

Insert

For the optimum loading of infusion flasks. Total area 84 x 84 mm.



- 20 compartments, dimension: 84 x 84 mm (bottle), 46 x 46 mm (neck)

EAN: 4002511115866 / Material number: 03808550 / Old Material Number: 69512901D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
PG 8535	•
PG 8536	•
PG 8562	•
Product category	
Insert for various bottles	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Number of compartments	20
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	113
External dimensions, net width in mm	445
External dimensions, net depth in mm	445
External dimensions, gross height in mm	18 15/16 (480)
External dimensions, gross width in mm	18 15/16 (480)
External dimensions, gross depth in in (mm)	5 1/8 (130)
Net weight in kg	4 (1.41)
Gross weight in kg	4 (1.78)