

E 134
Insert
for optimum loading of 210 slides.



- suitable for 210 slides
- Equipped with 26 compartments at 0 7/16 (11) in (mm) x 1 1/16 (26) in (mm), wire thickness 3 mm

EAN: 4002511115897 / Material number: 03808600 / Old Material Number: 69513401D

Product category	
Glassware Washers, laboratory	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PG 8582 CD	•
PG 8504	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 6011	•
PLW 6111	•
PLW 7111	•
PLW 8505	•
Product category	
Other inserts	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Number of compartments	26
Compartment width in mm	0 7/16 (11)
Compartment length in mm	1 1/16 (26)
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Slides [number]	210
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in in. (mm)	2 7/8 (73)
External dimensions, net width in in. (mm)	7 7/8 (200)
External dimensions, net depth in in. (mm)	17 9/16 (445)
External dimensions, gross height in in. (mm)	3 3/16 (80)
External dimensions, gross width in in. (mm)	9 1/16 (230)
External dimensions, gross depth in in. (mm)	18 1/8 (460)
Net weight in lb (kg)	4 (1.50)
Gross weight in lb (kg)	4 (1.61)