

E 135/1

Insert

for optimum loading of 19 baby bottles up to 110 ml.



- Suitable for 19 baby bottles x 110 ml
- Suitable for bottles sized 51 x 51 mm, neck: 45 x 45 mm
- For upper basket A 102 and lower basket A 151

EAN: 4002511115910 / Material number: 03808640 / Old Material Number: 69513502D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7826	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Insert for various bottles	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of baby bottles	•
Capacity	
Baby bottles [number]	19
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	5 3/8 (135) in (mm)
External dimensions, net width in mm	7 9/16 (192) in (mm)
External dimensions, net depth in mm	17 5/8 (447) in (mm)
External dimensions, gross height in mm	7 11/16 (195)
External dimensions, gross width in mm	8 1/8 (205)
External dimensions, gross depth in in (mm)	18 5/16 (465)
Net weight in kg	5 (2.02)
Gross weight in kg	5 (2.10)