

E 135/2

Insert

for optimum loading of 19 baby bottles up to 90 ml.

EAN: 4002511115927 / Material number: 03808650 / Old Material Number: 69513503D



- Suitable for 19 baby bottles x 90 ml
- Suitable for bottles sized 51 x 51 mm, neck: 42 x 42 mm
- For upper basket A 102 and lower basket A 151

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7826	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Insert for various bottles	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of baby bottles	•
Capacity	
Baby bottles [number]	19
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	4 15/16 (125) in (mm)
External dimensions, net width in mm	7 9/16 (192) in (mm)
External dimensions, net depth in mm	17 5/8 (447) in (mm)
External dimensions, gross height in mm	5 9/16 (140)
External dimensions, gross width in mm	9 5/16 (235)
External dimensions, gross depth in in (mm)	18 15/16 (480)
Net weight in kg	5 (1.97)
Gross weight in kg	5 (2.10)