

E 149

Insert

for optimum loading of 80 test tubes, 16 x 105 mm.



- Suitable for 80 test tubes up to 16 x 105 mm
- Additionally equipped with a lid
- 80 compartments at 0 3/4 (18) in (mm) x 0 3/4 (18) in (mm)
- Mesh width on the base 8 x 8 mm

EAN: 4002511116016 / Material number: 03808800 / Old Material Number: 69514901D

<b>Product category</b>	
Glassware Washers, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 6011	•
PLW 6111	•
PLW 8505	•
<b>Product category</b>	
Insert / module for test tubes	•
<b>Product characteristics</b>	
Material	Stainless steel
Colour	Stainless Steel
Number of compartments	80
Compartment width in mm	0 3/4 (18)
Compartment length in mm	0 3/4 (18)
<b>Application area</b>	
Reprocessing of test tubes	•
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Test tubes [number]	80
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	6 (152) in (mm)
External dimensions, net width in mm	7 7/8 (200) in (mm)
External dimensions, net depth in mm	7 7/8 (200) in (mm)
External dimensions, gross height in mm	6 3/4 (170)
External dimensions, gross width in mm	9 7/8 (250)
External dimensions, gross depth in in (mm)	9 7/8 (250)
Net weight in kg	3 (1.03)
Gross weight in kg	4 (1.59)