

E 331
Mobile unit
for the optimum loading of butyrometers.



- equipped with 39 injector nozzles

EAN: 4002511116115 / Material number: 03809350 / Old Material Number: 69533101D

Product category	
Glassware washers, laboratory	•
Machine affiliation	
PG 8535	•
PG 8536	•
Product category	
Other mobile units	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Reprocessing of butyrometers	•
Capacity	
Butyrometers [number]	39
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in in. (mm)	19 13/16 (502)
External dimensions, net width in in. (mm)	21 1/8 (535)
External dimensions, net depth in in. (mm)	20 3/8 (517)
External dimensions, gross height in in. (mm)	19 11/16 (500)
External dimensions, gross width in in. (mm)	21 5/16 (540)
External dimensions, gross depth in in. (mm)	20 7/8 (530)
Net weight in lb (kg)	12 (5.20)
Gross weight in lb (kg)	16 (7.10)
Accessory pack	
SD-B injector nozzle for butyrometer [number]	39
Rinse nozzle for powder dispenser [number]	1