

E 353

Clamps

Height-adjustable for optimum positioning of laboratory glassware on the nozzle.

EAN: 4002511117099 / Material number: 03809530 / Old Material Number: 69735301D



- Suitable for injector nozzle E 351
- Height adjustable locking
- Suitable for bottles, Erlenmeyer flasks and round bottomed flasks for example

Product category	
Glassware Washers, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
Product category	
Other accessories	•
Product characteristics	
Material	Wire, rilsan-coated, Stainless steel
Colour	Stainless Steel, White
Application area	
Reprocessing of laboratory glassware	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	3 15/16 (100) in (mm)
External dimensions, net width in mm	1 13/16 (45) in (mm)
External dimensions, net depth in mm	6 5/16 (160) in (mm)
External dimensions, gross height in mm	4 3/4 (120)
External dimensions, gross width in mm	3 3/16 (80)
External dimensions, gross depth in in (mm)	3 15/16 (100)
Net weight in kg	1 (0.03)
Gross weight in kg	1 (0.04)