

E 364

Insert

for the optimum loading of 36 wide-necked teats.



- Suitable for 36 wide-necked teats
- Equipped with 36 compartments at 1 3/4 (44) in (mm) x 1 3/4 (44) in (mm)
- Flip-up and lockable lid

EAN: 4002511116290 / Material number: 03809640 / Old Material Number: 69536401D

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7826	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Other inserts	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Number of compartments	36
Compartment width in mm	1 3/4 (44)
Compartment length in mm	1 3/4 (44)
Application area	
Reprocessing of baby bottles	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	77
External dimensions, net width in mm	215
External dimensions, net depth in mm	445
External dimensions, gross height in mm	3 3/16 (80)
External dimensions, gross width in mm	9 5/16 (235)
External dimensions, gross depth in in (mm)	18 1/8 (460)
Net weight in kg	4 (1.48)
Gross weight in kg	5 (1.91)