

E 374

Insert

for optimum loading of ENT instruments such as inhalation connections.



- 24 comp. approx. 45 x 45 mm, 27 compartments approx. 12 x 12 mm
- Mesh width on base: 1,7 mm
- For use in upper and lower basket

EAN: 4002511116412 / Material number: 03809800 / Old Material Number: 69537401D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
G 7831	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Insert for ENT instruments	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Number of compartments	51
Application area	
Reprocessing of ENT instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	2 1/2 (63)
External dimensions, net width in in. (mm)	6 13/16 (173)
External dimensions, net depth in in. (mm)	17 9/16 (445)
External dimensions, gross height in in. (mm)	3 (75)
External dimensions, gross width in in. (mm)	7 1/2 (190)
External dimensions, gross depth in in. (mm)	18 1/8 (460)
Net weight in lb (kg)	3 (1.20)
Gross weight in lb (kg)	3 (1.20)