

E 381

Mobile unit

for the optimum loading of intubation material.



- 20 nozzles Ø 4.0 x 30 mm, 5 nozzles Ø 2.5 x 30 mm (5 loose)
- 6 x E 466 injector nozzles for breathing bags Ø 8.0 x 333 mm
- 2 x E 431 injector nozzles for bellows, Ø 8.0 x 193 mm
- Equipped with connection for hot air drying
- Bracket for ML/2 for wash cart recognition / spray arm control

EAN: 4002511116474 / Material number: 03809890 / Old Material Number: 69538101D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
PG 8535	•
PG 8536	•
Product category	
Mobile unit for anaesthesia instruments	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing anaesthetic instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	19 13/16 (502)
External dimensions, net width in in. (mm)	21 1/8 (535)
External dimensions, net depth in in. (mm)	20 5/16 (515)
External dimensions, gross height in in. (mm)	21 5/16 (540)
External dimensions, gross width in in. (mm)	20 7/8 (530)
External dimensions, gross depth in in (mm)	19 11/16 (500)
Net weight in lb (kg)	15 (6.80)
Gross weight in lb (kg)	20 (8.70)