

E 408
 Mobile unit
 for the optimum loading of up to 96 pipettes.



- Suitable for 96 pipettes up to 450 mm
- Surface area 16 x 16 mm
- Complete with frame height 150 mm
- Equipped with connection for hot air drying

EAN: 4002511116689 / Material number: 03810130 / Old Material Number: 69540801D

Product category	
Glassware Washers, laboratory	•
Machine affiliation	
PG 8535	•
PG 8536	•
Product category	
Mobile unit with connection for hot-air drying	•
Mobile unit for pipettes	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Application area	
Reprocessing of pipettes / phials	•
Capacity	
Pipettes [number]	96
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	19 13/16 (502)
External dimensions, net width in in. (mm)	21 1/8 (535)
External dimensions, net depth in in. (mm)	20 3/8 (517)
External dimensions, gross height in in. (mm)	21 5/16 (540)
External dimensions, gross width in in. (mm)	20 7/8 (530)
External dimensions, gross depth in in. (mm)	19 11/16 (500)
Net weight in lb (kg)	13 (5.50)
Gross weight in lb (kg)	16 (7.10)