


E 413

Insert

for the optimum loading of aspirator hoses.



- For 6 System Sirona suction tubes
- For use in upper basket O177/1

EAN: 4002511124486 / Material number: 03810180 / Old Material Number: 68741301D

<b>Product category</b>	
Thermal disinfectors	•
<b>Machine affiliation</b>	
PG 8581	•
PG 8591	•
<b>Product category</b>	
Insert / module for dental instruments	•
<b>Product characteristics</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Application area</b>	
Reprocessing of dental instruments	•
<b>Capacity</b>	
Hoses [number]	6
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	205
External dimensions, net width in mm	390
External dimensions, net depth in mm	450
External dimensions, gross height in mm	9 7/8 (250)
External dimensions, gross width in mm	19 1/8 (485)
External dimensions, gross depth in in (mm)	19 1/2 (495)
Net weight in kg	4 (1.40)
Gross weight in kg	4 (1.41)