


A 13

Lid

For fixing delicate wash items.



- Lid in stainless steel
- 1 mm wire, 8 mm mesh size, 4 mm frames
- As replacement for E 103, E 104, E 105, and E 139 inserts

EAN: 4002511117204 / Material number: 03810200 / Old Material Number: 69745001D

Product category	
Glassware Washers, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 6011	•
PLW 6111	•
PLW 8505	•
Product category	
Lid	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in in. (mm)	0 13/16 (20)
External dimensions, net width in in. (mm)	7 7/8 (200)
External dimensions, net depth in in. (mm)	12 5/8 (320)
External dimensions, gross height in in. (mm)	0 13/16 (20)
External dimensions, gross width in in. (mm)	8 11/16 (220)
External dimensions, gross depth in in. (mm)	11 13/16 (300)
Net weight in lb (kg)	1 (0.14)
Gross weight in lb (kg)	1 (0.25)