


UTS

Utensil box

for the ideal cleaning of small parts.



- Equipped with lid
- 3 compartments with a size of ?? x ??

EAN: 4002511117143 / Material number: 03810260 / Old Material Number: 69741001D

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
G 7831	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Other inserts	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Number of compartments	3
Application area	
Reprocessing of MIS instruments	•
Reprocessing of surgical instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	93
External dimensions, net width in mm	102
External dimensions, net depth in mm	350
External dimensions, gross height in mm	4 3/8 (110)
External dimensions, gross width in mm	15 3/4 (400)
External dimensions, gross depth in in (mm)	23 7/8 (605)
Net weight in kg	2 (0.75)
Gross weight in kg	3 (1.00)