


ID 160

Injector nozzle

with plastic supports, Ø 0 3/16 (4) in (mm), length 160 mm.



- Suitable for standard mobile injector units and injector baskets
- Suitable for measuring flasks and test tubes, for example

EAN: 4002511116993 / Material number: 03810350 / Old Material Number: 69716001D

Product category	
Glassware Washers, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Injector nozzle	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Application area	
Reprocessing of laboratory glassware	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	2 3/16 (55)
External dimensions, net width in in. (mm)	2 3/16 (55)
External dimensions, net depth in in. (mm)	6 5/16 (160)
External dimensions, gross height in in. (mm)	2 3/16 (55)
External dimensions, gross width in in. (mm)	2 3/16 (55)
External dimensions, gross depth in in. (mm)	6 5/16 (160)
Diameter in in (mm)	0 3/16 (4)
Net weight in lb (kg)	1 (0.02)
Gross weight in lb (kg)	1 (0.02)