

E 118

Insert

for the optimum loading of half petri dishes.



- Suitable for 38 half plate petri dishes with Ø 100 mm
- equipped with 38 holders, height 70 mm, distance approx. 26 mm
- For use in upper basket

EAN: 4002511115767 / Material number: 03830270 / Old Material Number: 69511801D

<b>Product classification</b>	
Glassware Washers, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 6011	•
PLW 6111	•
PLW 8505	•
<b>Product category</b>	
Insert / module for half petri dishes	•
<b>Product characteristics</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Half petri dishes [number]	38
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	120
External dimensions, net width in mm	460
External dimensions, net depth in mm	445
External dimensions, gross height in mm	6 3/8 (160)
External dimensions, gross width in mm	19 (480)
External dimensions, gross depth in in (mm)	19 (480)
Net weight in kg	4 (1.54)
Gross weight in kg	5 (1.90)