

E 137

Insert

for optimum loading of 56 half petri dishes.



- Suitable for 56 half petri dishes with Ø 100 mm
- Equipped with 56 holders, at 2 13/16 (70) in (mm) height, spacing approx. 1 1/16 (26) in (mm)
- For use in lower basket
- For use with insert E 136

EAN: 4002511115958 / Material number: 03830290 / Old Material Number: 69513701D

Product category	
Glassware Washers, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
SlimLine laboratory glassware washers	•
Machines for cleaning industrial parts	•
Large-chamber laboratory glasswashers for cleaning industrial parts	•
SlimLine laboratory glasswashers for cleaning industrial parts	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PLW 8636 [IPC]	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 6011	•
PLW 6111	•
PLW 7111	•
PLW 7111 [IPC]	•
PLW 8505	•
Product category	
Insert / module for half petri dishes	•
Application area	
Reprocessing of laboratory glassware	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Distance between rails in mm	1 1/16 (26)
Rail height in mm	2 13/16 (70)
Capacity	
Half petri dishes [number]	56
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0

E 137

Insert

for optimum loading of 56 half petri dishes.

EAN: 4002511115958 / Material number: 03830290 / Old Material Number: 69513701D



Dimensions and weight	
External dimensions, net height in in. (mm)	19 1/8 (485)
External dimensions, net width in in. (mm)	17 9/16 (445)
External dimensions, net depth in in. (mm)	3 3/4 (95)
External dimensions, gross height in in. (mm)	5 9/16 (140)
External dimensions, gross width in in. (mm)	20 7/8 (530)
External dimensions, gross depth in in. (mm)	22 7/8 (580)
Net weight in lb (kg)	4 (2.03)
Gross weight in lb (kg)	8 (3.45)