

E 142

Insert

For holding different instruments.



- 1 mm wire, 5 mm mesh size, 5 mm frames
- equipped with 2 hinged handles
- Maximum load 10 kg

EAN: 4002511115972 / Material number: 03830300 / Old Material Number: 69514201D

Product category	
Glassware Washers, laboratory	•
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7824	•
G 7826	•
PG 8527	•
PG 8528	•
G 7831	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8583	•
PG 8583 CD	•
PG 8591	•
PG 8592	•
PG 8593	•
Product category	
Mesh tray	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Reprocessing of GYN instruments	•
Reprocessing of ENT instruments	•
Reprocessing of surgical instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	2 3/16 (55) in (mm)
External dimensions, net width in mm	10 1/16 (255) in (mm)
External dimensions, net depth in mm	18 15/16 (480) in (mm)
External dimensions, gross height in mm	3 (75)
External dimensions, gross width in mm	11 1/16 (280)
External dimensions, gross depth in in (mm)	20 5/16 (515)
Net weight in kg	2 (0.70)
Gross weight in kg	2 (0.86)