

AK 12
Insert
for holding different utensils.



- suitable for breast milk pumps, beakers & funnels, for example

EAN: 4002511115620 / Material number: 03830510 / Old Material Number: 69501201D

Product category	
Glassware Washers, laboratory	•
Washer-disinfectors, medical	•
Machines for cleaning industrial parts	•
Machine affiliation	
PG 8535	•
PG 8536	•
PLW 8636 [IPC]	•
PLW 8636 [LAB / LAB MON]	•
PG 8562	•
Product category	
Insert / module for beakers	•
Insert / module for wide-necked glassware	•
Other inserts	•
Application area	
Reprocessing of laboratory glassware	•
Reprocessing of different utensils	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in in. (mm)	4 13/16 (122)
External dimensions, net width in in. (mm)	8 7/8 (224)
External dimensions, net depth in in. (mm)	17 1/8 (434)
External dimensions, gross height in in. (mm)	4 7/8 (123)
External dimensions, gross width in in. (mm)	12 7/16 (315)
External dimensions, gross depth in in. (mm)	19 15/16 (505)
Net weight in lb (kg)	2 (0.74)
Gross weight in lb (kg)	2 (1.06)