


O 187

Upper basket

for the optimum loading of narrow necked glassware.



- Equipped with 34 injector nozzles
- Vertical clearance 6 3/4 (170) in (mm) mm\*
- Powder dispensing not possible

\* Depending on the laboratory glassware the max. vertical clearance can be lower.

EAN: 4002511808881 / Material number: 03881640 / Old Material Number: 69118701D

<b>Product category</b>	
Glassware Washers, laboratory	•
<b>Machine affiliation</b>	
PG 8535	•
PG 8536	•
<b>Product category</b>	
Injector basket	•
Upper basket	•
<b>Product characteristics</b>	
Material	Stainless steel
Colour	Stainless Steel
Maximum vertical clearance in mm	6 3/4 (170)
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Narrow necked glassware [number]	34
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in in. (mm)	9 7/8 (250)
External dimensions, net width in in. (mm)	18 3/4 (475)
External dimensions, net depth in in. (mm)	20 15/16 (531)
External dimensions, gross height in in. (mm)	10 11/16 (270)
External dimensions, gross width in in. (mm)	20 11/16 (525)
External dimensions, gross depth in in (mm)	21 5/16 (540)
Net weight in lb (kg)	13 (5.60)
Gross weight in lb (kg)	15 (6.50)
<b>Accessory pack</b>	
E 351 injector nozzle, Ø 4 mm, length 160 mm [number]	34
E 353 clip for injector nozzle E 351 [number]	34