

E 416

Insert

for optimum loading of 6 single and two-piece specula.



- Suitable for 6 single and two-part specula
- 7 holders 40 mm apart
- For use in upper (only with single-part specula) and lower basket

EAN: 4002511803824 / Material number: 03888650 / Old Material Number: 69541601D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
G 7831	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
PWD 8534	•
Product category	
Insert for gynaecology instruments	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of GYN instruments	•
Capacity	
GYN specula [number]	6
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	6 3/16 (157)
External dimensions, net width in in. (mm)	7 1/16 (178)
External dimensions, net depth in in. (mm)	11 (279)
External dimensions, gross height in in. (mm)	8 5/16 (210)
External dimensions, gross width in in. (mm)	7 7/8 (200)
External dimensions, gross depth in in (mm)	14 3/8 (365)
Net weight in lb (kg)	2 (0.69)
Gross weight in lb (kg)	3 (0.97)