

MT
Trolleys
for loading & unloading washer-&disinfectors.

EAN: 4002511820005 / Material number: 03893140 / Old Material Number: 69200102D



- For the storage and transport of baskets and inserts
- Divided into 4 levels that can be adjusted by 102.5 mm
- Moveable with 4 lockable castors
- Loading volume 549 x 599 mm

Product category	
Glassware Washers, laboratory	•
Washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
G 7831	•
PG 8535	•
PG 8536	•
PG 8582	•
PG 8582 CD	•
PG 8583	•
PG 8583 CD	•
PG 8592	•
PG 8593	•
Product category	
Storage / transport system	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	78 3/16 (1985)
External dimensions, net width in in. (mm)	26 1/8 (662)
External dimensions, net depth in in. (mm)	24 5/16 (616)
External dimensions, gross height in in. (mm)	79 15/16 (2030)
External dimensions, gross width in in. (mm)	26 13/16 (680)
External dimensions, gross depth in in. (mm)	26 13/16 (680)
Net weight in lb (kg)	51 (23.00)
Gross weight in lb (kg)	64 (29.00)