


A 14

Lid

For fixing delicate wash items.



- Lid in stainless steel
- 7 x 7 mm perforations, 3 mm slat
- For AK 12 insert

EAN: 4002511860995 / Material number: 03981970 / Old Material Number: 69745002D

Product category	
Glassware washers, laboratory	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8535	•
PG 8536	•
PG 8562	•
Product category	
Lid	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Reprocessing of different utensils	•
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in in. (mm)	0 13/16 (20)
External dimensions, net width in in. (mm)	8 5/16 (210)
External dimensions, net depth in in. (mm)	8 5/16 (210)
External dimensions, gross height in in. (mm)	1 3/16 (30)
External dimensions, gross width in in. (mm)	9 1/2 (240)
External dimensions, gross depth in in. (mm)	9 7/8 (250)
Net weight in lb (kg)	1 (0.20)
Gross weight in lb (kg)	1 (0.29)