


A 16
Lid
For fixing delicate wash items.



- 1 mm wire, 5 mm mesh size, 4 mm frames

EAN: 4002511860971 / Material number: 03981990 / Old Material Number: 69745004D

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
G 7831	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Lid	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	0 3/4 (19) in (mm)
External dimensions, net width in mm	9 1/16 (230) in (mm)
External dimensions, net depth in mm	8 1/2 (215) in (mm)
External dimensions, gross height in mm	1 5/8 (40)
External dimensions, gross width in mm	9 7/8 (250)
External dimensions, gross depth in in (mm)	10 1/16 (255)
Net weight in kg	1 (0.15)
Gross weight in kg	1 (0.30)