

E 442

Irrigation sleeves

for the optimum cleaning of MIS instruments.

EAN: 4002511908208 / Material number: 04203830 / Old Material Number: 69744201D



- Suitable for instruments with \varnothing 4 – 8 mm, screw connection
- Length 121 mm mm, \varnothing 0 1/4 (6) in (mm) mm
- Cap, opening \varnothing 6 mm
- Mounting clip (mat. no. 4174850)

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
Product category	
Irrigation sleeve	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Application area	
Reprocessing of MIS instruments	•
Reprocessing of lumened instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	121
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	0 1/2 (12)
External dimensions, gross width in mm	3 1/4 (80)
External dimensions, gross depth in in (mm)	4 3/4 (120)
Diameter in in (mm)	0 1/4 (6)
Net weight in kg	1 (0.03)
Gross weight in kg	1 (0.03)