

E 448

Silicone hose

for the optimum cleaning of various instruments.



- Luer lock adapter, male
- Length 11 13/16 (300) in (mm)
mm mm, Ø 0 7/16 (10) in (mm) mm
- Jets, thread 8 x 1 mm

EAN: 4002511908475 / Material number: 04203880 / Old Material Number: 69744801D

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Other accessories	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Transparent
Application area	
Reprocessing of MIS instruments	•
Reprocessing of ophthalmology instruments	•
Reprocessing of lumened instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	11 13/16 (300) in (mm)
External dimensions, net width in mm	0 3/8 (8) in (mm)
External dimensions, net depth in mm	0 3/8 (8) in (mm)
External dimensions, gross height in mm	0 7/16 (10)
External dimensions, gross width in mm	6 3/4 (170)
External dimensions, gross depth in in (mm)	4 15/16 (125)
Diameter in in (mm)	0 7/16 (10)
Net weight in kg	1 (0.02)
Gross weight in kg	1 (0.02)