


E 456
Spring stay
for the optimum hold of open MIS instruments.



EAN: 4002511908536 / Material number: 04203960 / Old Material Number: 69745601D

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Holder	•
Other accessories	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of MIS instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	3 (76) in (mm)
External dimensions, net width in mm	0 1/8 (2) in (mm)
External dimensions, net depth in mm	5 11/16 (144) in (mm)
External dimensions, gross height in mm	0 1/2 (12)
External dimensions, gross width in mm	3 3/16 (80)
External dimensions, gross depth in in (mm)	4 3/4 (120)
Net weight in kg	1 (0.00)
Gross weight in kg	1 (0.01)