

E 458

Insert

for the optimum loading of screw-cap teats.



- Suitable for 36 screw connection teats
- 36 compartments at ?? x ??
- equipped with a flip-up and lockable lid

EAN: 4002511919358 / Material number: 04235200 / Old Material Number: 69545801D

<b>Product category</b>	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7823	•
G 7824	•
G 7826	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
<b>Product category</b>	
Other inserts	•
<b>Product characteristics</b>	
Material	Stainless steel
Colour	Stainless Steel
Number of compartments	36
<b>Application area</b>	
Reprocessing of baby bottles	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in in. (mm)	2 1/2 (63)
External dimensions, net width in in. (mm)	8 1/2 (215)
External dimensions, net depth in in. (mm)	17 9/16 (445)
External dimensions, gross height in in. (mm)	17 3/4 (450)
External dimensions, gross width in in. (mm)	3 3/16 (80)
External dimensions, gross depth in in. (mm)	8 11/16 (220)
Net weight in lb (kg)	4 (1.50)
Gross weight in lb (kg)	4 (1.60)