

E 471

Spring clip

for injector jets with a Ø of 2.5 mm.



EAN: 4002512008464 / Material number: 04451930 / Old Material Number: 69747101D

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Other accessories	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing anaesthetic instruments	•
Reprocessing of MIS instruments	•
Reprocessing of ENT instruments	•
Reprocessing of lumened instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	0 1/2 (12)
External dimensions, gross width in mm	3 3/16 (80)
External dimensions, gross depth in in (mm)	4 3/4 (120)
Net weight in kg	1 (0.00)
Gross weight in kg	1 (0.00)