

E 800
Insert
for optimum loading of 3 mesh or kidney trays.



- Equipped with 4 holders
- For use in upper and lower basket

EAN: 4002512087773 / Material number: 04637720 / Old Material Number: 69580001D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
G 7831	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Insert for mesh trays	•
Insert for kidney dishes / basins	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of kidney dishes	•
Capacity	
Kidney dishes [number]	3
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	165
External dimensions, net width in mm	140
External dimensions, net depth in mm	290
External dimensions, gross height in mm	7 1/2 (190)
External dimensions, gross width in mm	7 11/16 (195)
External dimensions, gross depth in in (mm)	16 15/16 (430)
Net weight in kg	1 (0.29)
Gross weight in kg	2 (0.49)