

E 941
Mobile unit
for modules on 2 levels.



- Suitable for up to 2 modules per level depending on the size of the load
- Coupling of the mobile unit to the water and drying air supply
- The clearances for both levels are different

EAN: 4002512273411 / Material number: 04812530 / Old Material Number: 69594101D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
PG 8527	•
PG 8528	•
Product category	
Mobile unit for orthopaedic instruments	•
Other mobile units	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	16 5/8 (421)
External dimensions, net width in in. (mm)	24 3/8 (619)
External dimensions, net depth in in. (mm)	31 1/8 (790)
External dimensions, gross height in in. (mm)	27 9/16 (700)
External dimensions, gross width in in. (mm)	31 1/2 (800)
External dimensions, gross depth in in (mm)	47 1/4 (1200)
Net weight in lb (kg)	36 (16.30)
Gross weight in lb (kg)	58 (26.00)
Accessory pack	
Magnetic strip for automatic program recognition (incl. magnets) [number]	1