

DOS S

Dispenser pump rack

With slanted base for the complete emptying of the liquid media containers.



- Suitable for the installation of the dispenser pump casing DOS G2
- Suitable for the installation of the dispenser pump unit DOS GP3
- With four swivel castors for maximum mobility

EAN: 4002512161107 / Material number: 04814240 / Old Material Number: 59170031D

Product classification	
Washing machines	•
Machine / accessory category	
PW 413	•
PW 418	•
PW 5064	•
PW 5062	•
PW 5055	•
PW 5065	•
PW 5064 MOP STAR 60	•
PW 5070	•
PW 5082	•
PW 5084 MOP STAR 80	•
PW 5104 MOP STAR 100	•
PW 5105 Vario	•
PW 5134 MOP STAR 130	•
PW 6055	•
PW 6065	•
PW 6055 Vario	•
PW 6065 Vario	•
PW 6065 Plus	•
PW 6068 Plus	•
PW 6080	•
PW 6088	•
PW 6107	•
PW 6137	•
PW 5136	•
PW 6163	•
PW 6167	•
PW 6207	•
PW 6241	•
PW 6243	•
PW 6321	•
PW 6323	•
PW 811	•
PW 814	•
PW 818	•
PW 9055	•
PW 9065	•
PWT 6089	•
Product category	
Dispenser system	•
Product characteristics	
Electrical connection	OHNE
Material	Powder coated galvanised steel
Colour	Iron grey
Dimensions and weight	
External dimensions, net height in mm	40 5/8 (1030)
External dimensions, net width in mm	39 3/8 (1000)

DOS S

Dispenser pump rack

With slanted base for the complete emptying of the liquid media containers.

EAN: 4002512161107 / Material number: 04814240 / Old Material Number: 59170031D



External dimensions, net depth in mm	14 5/8 (371)
External dimensions, gross height in mm	40 5/8 (1030)
External dimensions, gross width in mm	51 1/4 (1300)
External dimensions, gross depth in mm	16 3/8 (413)
Net weight [kg]	20 (9.00)
Gross weight in kg	23 (10.00)