

E 466
 Injector nozzle
 for the optimum loading of breathing bags 8 x 333 mm.



EAN: 4002512223744 / Material number: 04966620 / Old Material Number: 69746601D

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Other accessories	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing anaesthetic instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	0 (0)
External dimensions, net width in in. (mm)	0 (0)
External dimensions, net depth in in. (mm)	13 1/8 (333)
External dimensions, gross height in in. (mm)	5 9/16 (140)
External dimensions, gross width in in. (mm)	10 1/4 (260)
External dimensions, gross depth in in. (mm)	5 15/16 (150)
Diameter in in (mm)	0 3/8 (8)
Net weight in lb (kg)	2 (0.85)
Gross weight in lb (kg)	3 (1.00)