

E 930

Insert

for optimum loading of up to 36 theatre shoes up to size 44.



- For use in lower level of mobile unit E 975/2
- Equipped with 36 holders, height 280 mm

EAN: 4002512257145 / Material number: 05001360 / Old Material Number: 69593001D

Product category	
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
PG 8527	•
PG 8528	•
Product category	
Insert for OR shoes	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of OR shoes	•
Capacity	
OR shoes [number]	36
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	289
External dimensions, net width in mm	586
External dimensions, net depth in mm	752
External dimensions, gross height in mm	13 7/8 (350)
External dimensions, gross width in mm	30 3/8 (770)
External dimensions, gross depth in in (mm)	39 3/8 (1000)
Net weight in kg	9 (3.65)
Gross weight in kg	13 (5.70)