

E 484
Insert
for holding different utensils



- Wire mesh 1.4 mm, mesh width 8 mm
- Can be equipped with holders: E 485, E 486, E 487, E 488, E 489

EAN: 4002512398718 / Material number: 05109190 / Old Material Number: 69548401D

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Insert for kidney dishes / basins	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of OR shoes	•
Reprocessing of basins	•
Reprocessing of kidney dishes	•
Reprocessing of surgical instruments	•
Reprocessing of different utensils	•
Capacity	
OR shoes [number]	16
Basins [number]	4
Kidney dishes [number]	9
Breathing masks [number]	9
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	5 15/16 (150)
External dimensions, net width in mm	18 9/16 (470)
External dimensions, net depth in mm	18 15/16 (480)
External dimensions, gross height in mm	22 1/16 (560)
External dimensions, gross width in mm	7 1/8 (180)
External dimensions, gross depth in in (mm)	21 5/16 (540)
Net weight in kg	7 (2.97)
Gross weight in kg	9 (3.78)