

E 486
Bracket
for the optimum loading of basins.



- Suitable for 1 basins

EAN: 4002512398732 / Material number: 05109230 / Old Material Number: 69548601D

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Insert for kidney dishes / basins	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of basins	•
Capacity	
Basins [number]	1
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	275
External dimensions, net width in mm	464
External dimensions, net depth in mm	10
External dimensions, gross height in mm	18 3/4 (475)
External dimensions, gross width in mm	18 3/8 (465)
External dimensions, gross depth in in (mm)	1 (25)
Net weight in kg	1 (0.21)
Gross weight in kg	1 (0.43)