

E 487

Bracket

for optimum loading of 4 theatre shoes up to a height of 280 mm.

EAN: 4002512398749 / Material number: 05109250 / Old Material Number: 69548701D



Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Insert for OR shoes	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of OR shoes	•
Capacity	
OR shoes [number]	4
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	280
External dimensions, net width in mm	464
External dimensions, net depth in mm	10
External dimensions, gross height in mm	18 5/8 (470)
External dimensions, gross width in mm	18 5/8 (470)
External dimensions, gross depth in in (mm)	1 (25)
Net weight in kg	1 (0.37)
Gross weight in kg	2 (0.52)