


E 488
Bracket
for the optimum loading of respiratory masks.



- Suitable for 3 respirator masks

EAN: 4002512398756 / Material number: 05109260 / Old Material Number: 69548801D

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Other inserts	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of different utensils	•
Capacity	
Breathing masks [number]	3
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	1 1/4 (31) in (mm)
External dimensions, net width in mm	9 1/2 (240) in (mm)
External dimensions, net depth in mm	8 5/16 (210) in (mm)
External dimensions, gross height in mm	3 (75)
External dimensions, gross width in mm	11 1/16 (280)
External dimensions, gross depth in in (mm)	20 1/8 (510)
Net weight in kg	1 (0.27)
Gross weight in kg	1 (0.39)