

E 713

Mobile unit

for loading of 2 DIN containers including lids as well as 4 mesh trays on 3 levels.



- Equipped with 2 integrated spray arms
- The loading areas for the individual levels are different
- Magnetic strip for mobile unit recognition (without magnets ML)
- Equipped with a connection for 5 hollow instruments
- Equipped with connection for hot air drying

EAN: 4002512401944 / Material number: 05221550 / Old Material Number: 69571301D

Product category	
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7825	•
G 7826	•
Product category	
Mobile unit for inserts	•
Mobile unit for containers	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of surgical instruments	•
Reprocessing of containers	•
Capacity	
DIN mesh trays [number]	4
SPRI mesh trays [number]	2
DIN containers, height 15 cm [number]	2
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	20 7/8 (530)
External dimensions, net width in in. (mm)	22 11/16 (575)
External dimensions, net depth in in. (mm)	23 5/8 (600)
External dimensions, gross height in in. (mm)	22 1/16 (560)
External dimensions, gross width in in. (mm)	25 1/4 (640)
External dimensions, gross depth in in (mm)	26 13/16 (680)
Net weight in lb (kg)	22 (9.70)
Gross weight in lb (kg)	23 (10.00)
Accessory pack	
Magnetic strip for automatic program recognition (without magnets) [number]	1