


E 475
Holding bar
for stabilising sleeves.



- Suitable for MIS inserts

EAN: 4002512401234 / Material number: 05226040 / Old Material Number: 69747501D

Product classification	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Holder	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of MIS instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	110
External dimensions, net width in mm	155
External dimensions, net depth in mm	500
External dimensions, gross height in mm	4 3/8 (110)
External dimensions, gross width in mm	6 5/16 (160)
External dimensions, gross depth in in (mm)	20 1/8 (510)
Net weight in kg	1 (0.10)
Gross weight in kg	1 (0.22)