

E 496

Injector nozzle

with retaining clamp for intubation material.



EAN: 4002512365970 / Material number: 05235490 / Old Material Number: 69749601D

<b>Product category</b>	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
<b>Product category</b>	
Injector nozzle	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
<b>Application area</b>	
Reprocessing anaesthetic instruments	•
Reprocessing of MIS instruments	•
Reprocessing of lumened instruments	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in in. (mm)	4 3/4 (120)
External dimensions, net width in in. (mm)	0 (0)
External dimensions, net depth in in. (mm)	0 (0)
External dimensions, gross height in in. (mm)	4 3/4 (120)
External dimensions, gross width in in. (mm)	1 3/16 (30)
External dimensions, gross depth in in. (mm)	1 3/16 (30)
Diameter in in (mm)	0 3/16 (4)
Net weight in lb (kg)	2 (0.75)
Gross weight in lb (kg)	2 (0.90)