

E 730  
 Insert  
 for the optimum loading of 20 theatre shoes.



EAN: 4002512423120 / Material number: 05271020 / Old Material Number: 69573001D

<b>Product category</b>	
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7825	•
G 7826	•
<b>Product category</b>	
Insert for OR shoes	•
<b>Product characteristics</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Application area</b>	
Reprocessing of OR shoes	•
<b>Capacity</b>	
OR shoes [number]	20
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in in. (mm)	11 7/16 (290)
External dimensions, net width in in. (mm)	18 9/16 (470)
External dimensions, net depth in in. (mm)	21 5/16 (540)
External dimensions, gross height in in. (mm)	16 9/16 (420)
External dimensions, gross width in in. (mm)	21 5/16 (540)
External dimensions, gross depth in in (mm)	23 5/8 (600)
Net weight in lb (kg)	20 (8.70)
Gross weight in lb (kg)	23 (10.00)