


E 984
Insert
for a wide variety of utensils.



- Can be equipped with different holders

EAN: 4002512431323 / Material number: 05337750 / Old Material Number: 69598401D

Product category	
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
PG 8527	•
PG 8528	•
Product category	
Insert for kidney dishes / basins	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of OR shoes	•
Reprocessing of basins	•
Reprocessing of kidney dishes	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	371
External dimensions, net width in mm	150
External dimensions, net depth in mm	572
External dimensions, gross height in mm	15 3/8 (390)
External dimensions, gross width in mm	6 3/8 (160)
External dimensions, gross depth in in (mm)	23 1/4 (590)
Net weight in kg	7 (2.95)
Gross weight in kg	8 (3.30)