

E 985
Bracket
for optimum loading of 2 kidney trays.



- For use in insert E 984

EAN: 4002512431330 / Material number: 05337760 / Old Material Number: 69598501D

Product category	
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
PG 8528	•
Product category	
Insert for kidney dishes / basins	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of kidney dishes	•
Capacity	
Kidney dishes [number]	2
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	2 7/8 (70)
External dimensions, gross width in mm	11 1/8 (280)
External dimensions, gross depth in in (mm)	20 (505)
Net weight in kg	1 (0.23)
Gross weight in kg	1 (0.43)