

E 986
Bracket
for the optimum loading of 5 basins.



- For use in insert E 984

EAN: 4002512431347 / Material number: 05337770 / Old Material Number: 69598601D

| | |
|---|-----------------|
| Product category | |
| Large-capacity washer-disinfectors, medical | • |
| Machine affiliation | |
| PG 8527 | • |
| PG 8528 | • |
| Product category | |
| Insert for kidney dishes / basins | • |
| Product characteristics | |
| Material | Stainless steel |
| Colour | Stainless Steel |
| Application area | |
| Reprocessing of basins | • |
| Capacity | |
| Basins [number] | 1 |
| Standard electrical connection | |
| Number of phases | 0 |
| Voltage in V | 0 |
| Frequency in Hz | 0 |
| Total rated load in kW | 0.00 |
| Fuse rating in A | 0 |
| Dimensions and weight | |
| External dimensions, net height in mm | 0 |
| External dimensions, net width in mm | 0 |
| External dimensions, net depth in mm | 0 |
| External dimensions, gross height in mm | 1 (25) |
| External dimensions, gross width in mm | 13 (330) |
| External dimensions, gross depth in in (mm) | 17 1/4 (435) |
| Net weight in kg | 1 (0.19) |
| Gross weight in kg | 1 (0.36) |