

E 987
Bracket
for optimum loading of 15 theatre shoes up to size 44.



- For use in insert E 984

EAN: 4002512431354 / Material number: 05337800 / Old Material Number: 69598701D

Product category	
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
PG 8527	•
PG 8528	•
Product category	
Insert for OR shoes	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of OR shoes	•
Capacity	
OR shoes [number]	15
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	0 (0)
External dimensions, net width in in. (mm)	0 (0)
External dimensions, net depth in in. (mm)	0 (0)
External dimensions, gross height in in. (mm)	1 3/16 (30)
External dimensions, gross width in in. (mm)	18 1/8 (460)
External dimensions, gross depth in in (mm)	18 3/4 (475)
Net weight in lb (kg)	1 (0.23)
Gross weight in lb (kg)	1 (0.40)